




EPO Analysis on SRT[®] SEC-300



Column: SRT[®] SEC-300 (5 μ m, 300 \AA , 7.8 x 300 mm)
Mobile Phase: $\text{Na}_2\text{HPO}_4\text{-KH}_2\text{PO}_4\text{-NaCl}$
Flow Rate: 0.7 mL/min
Detection: UV 280 nm
Temperature: Ambient
Injection Volume: 20.0 μ L
Sample: EPO (Recombinant Human Erythropoietin) (MW 60-70KD)

Keywords: size exclusion, SRT, SRT 300, protein, Recombinant Human Erythropoietin